

Why Bodyguards Glove Range

Despite the introduction of hand-washing initiatives, the level of hand-washing cannot be guaranteed to reach 100% compliance. Gloves, traditional glove boxes and dispensers can themselves be a vehicle for contamination!

Bodyguards SafeDon™ Gloves are dispensed individually from the bottom of the box. This unique dispensing mechanism ensures that the user does not touch the gloves within the box or the box itself whilst donning gloves, this helps to minimise cross contamination.

Dispensing any glove type "Cuff 1st" helps to reduce the contamination normally caused by the traditional glove dispensing methods by which the user can touch the glove box with dirty hands or remove gloves from the box by the fingers or palm areas.

The Bodyguards SafeDon™ glove box is securely contained within a dispensing unit that has an anti-microbial protective additive which helps to prevent the build up and growth of organisms on the glove dispenser.

The Bodyguards SafeDon™ system provides an effective infection control solution for use where hygiene is essential.

With appropriate use of gloves and the adherence to hand hygiene policies the principle of gloves dispensed "Cuff 1st" from a dispenser could potentially reduce the risk of cross contamination.

Reduce the Risk of Cross Contamination by utilising our New Unique Hygienic Glove Dispensing System...

It is widely acknowledged that hands are the main contamination route of infection. The current way of donning a disposable glove means the fingers, palms and thumbs of gloves in a traditional glove box are touched to allow removal of a glove from the box.

The Bodyguards SafeDon™ concept has been launched to provide a quick, hygienic and economic glove dispensing system for medical environments and has been specially designed to help combat the problem of cross contamination of gloves, glove boxes and dispensers when donning disposable gloves.

KEY BENEFITS

An innovative glove dispensing system featuring many key benefits to the user including:

NEW CUFF 1ST DISPENSING SYSTEM

Hygienic system for dispensing and donning gloves. Each glove is dispensed individually and can be donned hygienically "Cuff 1st" without touching the fingers and palm areas of the glove, this

helps to reduce the potential risk of cross contamination.

TIME EFFICIENT

Fast, easy dispensing saves time and patients waiting.

COST EFFECTIVE

Minimise glove wastage as gloves are dispensed individually and disused gloves cannot be put back into the box which helps further reduce cross contamination.

ERGONOMICALLY DESIGNED DISPENSER

The Bodyguards SafeDon™ acrylic wall mounted dispenser keeps gloves away from contaminated surfaces and can be wiped down to maintain hygiene. The acrylic dispenser features an anti-microbial additive and has been designed without joins to avoid harbouring bacteria and trapping moisture.

FLEXIBLE GLOVE OPTIONS

Each compact cardboard dispensing box contains 100 nitrile, latex or synthetic gloves (X/large glove boxes contain 90 gloves except for Vinyl). When the box of gloves is finished, simply remove the empty box and replace with a new box.

“The future is in your hands, you can help reduce the spread of infection and the associated financial burden on resources”